ARTI PLASTER A

EXTRA ADHESIVE

Artikulationsgips der Klasse III

Indikation

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Articulating plaster class III

Indication

- High adhesive strength
- Very low expansion
- Short processing time

Processing

100 g powder : 30 ml distilled / demineralised water

The amount of water can be increased or decreased by up to 1,5 ml. With less liquid the plaster will be harder and the working time will be shorter. Using more liquid has the reverse effect.

Please make sure only to mix enough plaster required for two to three articulations (this is important for crystallization!).

Empty the plaster into the mixing beaker and add water. Stir slowly by hand until powder and liquid are well mixed. Mix with spatula for another 20 - 30 seconds.

Important to follow the indicated mixing time and mixing ratio:

The longer the pulp is mixed and the thicker its consistency, the lower its adhesive force.

Technical data: (applies for a mixing ratio of 100 g : 30 ml)

Brinell hardness after 2 hours (MPa; N/mm ²).....	80
Brinell hardness after 24 hours (MPa; N/mm ²).....	90
Compressive strength after 1 hour (MPa; N/mm ²).....	29 (4,205 psi)
Setting expansion (%).....	< 0,03
End of setting (min).....	5
Working time (min).....	2,0 - 2,5

Shelf life

At least 24 months, if closed and stored under dry conditions.

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Articulating plaster class III

Indication

- Very low expansion
- Short processing time

Processing

100g powder: 30 ml distilled / demineralised water

The amount of water can be increased or decreased by up to 1,5 ml. With less liquid the plaster will be harder and the working time will be shorter. Using more liquid has the reverse effect. Please make sure only to mix enough plaster required for two to three articulations (this is important for crystallization!). Empty the plaster into the mixing beaker and add water. Stir slowly by hand until powder and liquid are well mixed. Mix with spatula for another 20 - 30 seconds.

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